


Booklet creator for medical image distribution

Dicom Paper

Dicom Paper is a versatile paper report creator for high volume imaging centers that allows users to create simple printed documents or booklets from anywhere on the network for distribution to patients.

Supporting the entire personalization of the document, Dicom Paper reduces the cost of report compilation, improves operational workflow and provides a highly professional product for patients and long term management.

Dicom Paper receives DICOM studies from multiple DICOM sources (Modalities, Servers, Viewers or Workstations) and creates a document that is printed with a personalized cover and back page, personalized header and footer, report page and the selected images for the report.



Included features:

- Touch interface for easy management
- Send studies directly from Modality
- Send studies from DICOM Viewer
- Re-print previous printed study
- Personalize cover page, and back page
- Personalize header and footer
- Personalize document with patient information included in DICOM header information
- Create profiles for printing documents from certain Modalities or Calling AE_Titles
- Highly customizable, flexible and scalable
- Settings blocked for administrators

Features for partners:

- Personalized Application with your own brand
- Include logo of brand and company
- Second-line support
- Online remote support
- Online private area:
 - Request licenses
 - Create new client accounts
 - Activate demo applications
 - Cancel and suspend applications

Grow you business. Become a Dicom Server Partner

www.dicomservers.com - info@seamlink.com - Seamlink, lda