



## Sincronization server for digital medical images

### Dicom Sync

Dicom Sync server is a versatile tool for high volume imaging centers that allows IT manager to setup rules for DICOM image sincronization between multiple equipments.

With a highly robust performance, Dicom Sync server reduces the cost of managing the flow of medical images inside and between imaging centers. It also provides tools for imaging backup.

Dicom Sync server receives DICOM studies from multiple DICOM sources (Modalities, Servers, Viewers or Workstations) and rediretcs them to a destination or stores them for backup. The server has the ability to establish rules so that it knows the destination of each type of image.

### Our services:

- Remote Diagnostics
- Onsite Support
- Installation
- Training
- Software Upgrades
- Disaster Recovery
- 1° year warranty

### Included features:

- Touch interface for easy management
- Send studies directly from Modality
- Send studies from DICOM Viewer
- Do Query & Retreive to PACS server
- Create rules to rewrite DICOM Datasets
- Create rules to determine the destination of each kind of study
- Highly customizable, flexible and scalable
- Settings blocked for administrators

### Features for partners:

- Personalized Application with your own brand
- Include logo of brand and company
- Second-line support
- Online remote support
- Online private area:
  - Request licenses
  - Create new client accounts
  - Activate demo applications
  - Cancel and suspend applications

Grow you business. Become a Dicom Server Partner

[www.dicomservers.com](http://www.dicomservers.com) - [info@seamlink.com](mailto:info@seamlink.com) - Seamlink, lda